




PAUL HEWITT  
**PHREP MEASURING TAPE**

Fits for all PAUL HEWITT PHREPS.

Make finding your fitting wristband size as easy as 1-2-3. With the PHREP measuring tape, we help you find the perfect size of your PAUL HEWITT PHREP in just four steps.



**1. Print it out**

Print as DIN A4 in original size.  
(actual size or 100%)



**2. Cut it out**

Cut out the measuring tape along the dotted line.



**3. Take the measurement**

Fit it around your wrist and read off the matching size on the arrows. If you are pending between two sizes please select the next larger PHREP.



**4. Order**

Order your desired PHREP at  
[www.paul-hewitt.com](http://www.paul-hewitt.com).

Size table (applicable for all PAUL HEWITT PHREPS)

XS	S	M	L	XL	XXL	3XL
16,0 cm/ 6.3"	17,0 cm/ 6.7"	18,0 cm/ 7.1"	19,0 cm/ 7.5"	20,0 cm/ 7.9"	21,0 cm/ 8.3"	22,0 cm/ 8.7"

**IMPORTANT**

Check, if the printout is **on the right scale**. To do so, simply check the measurement using a ruler, to see if the line for the scale is exactly 21,0 cm/8.3" long. Please **add 1 cm to the wrist circumference** if you use a standard measure tape, as the stated length in cm is for the whole PHREP.



**Reference size**

To verify the accuracy of the document printout scale please place a 1 Euro coin on the contour model.

